



GABBR2 (Phospho Ser893) rabbit pAb

Catalog No	YP-Ab-16367
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB; ELISA
Gene Name	GABBR2 GPR51 GPRC3B
Immunogen	Synthesized peptide derived from human GABBR2 (Phospho Ser893)
Specificity	This antibody detects endogenous levels of Human,Mouse,Rat GABBR2 (Phospho Ser893)
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	BcDNA:GH7312;CG676;CT2836;D Gaba2;FLJ36928;G protein coupled receptor 51;G-protein coupled receptor 51;GAB B R2;GABA B R2;GABA B receptor 2;GABA-B receptor 2;GABA-B-R2;GABA-BR2;GABAB R2;GABABR 2;GABABR2;GABB R2;GABBR 2;Gabbr2;GABR2_HUMAN;Gamma aminobutyric acid B receptor 2;Gamma aminobutyric acid GABA B receptor 2;Gamma aminobutyric acid type B receptor subunit 2;Gamma-aminobutyric acid type B receptor subunit 2;Gb 2;Gb2;GH7312;GPR 51;GPR51;GPRC 3B;GPRC3B;HG 2;HG2;HRIHFB299;Metabotropic GABA B receptor subtype 2;OTTHUMP21776;OTTHUMP63797;R2 SUBUNIT;GABBR2 GPR51 GPRC3B;GABBR2 (Phospho Ser893);Gamma-aminobutyric acid type B receptor subunit 2 (GABA-B receptor 2;HG20)
Observed Band	105kD
Calculated Molecular Weight	105kD
Cell Pathway	Cell membrane ; Multi-pass membrane protein . Cell junction, synapse, postsynaptic cell membrane ; Multi-pass membrane protein . Coexpression of GABBR1 and GABBR2 is required for GABBR1 maturation and transport to the plasma membrane. In contrast, GABBR2 does not depend on GABBR1 for transport to the cell membrane. .



Tissue Specificity	Highly expressed in brain, especially in cerebral cortex, thalamus, hippocampus, frontal, occipital and temporal lobe, occipital pole and cerebellum, followed by corpus callosum, caudate nucleus, spinal cord, amygdala and medulla (PubMed:10087195, PubMed:10328880, PubMed:10727622, PubMed:9872744) . Weakly expressed in heart, testis and skeletal muscle (PubMed:10087195, PubMed:10727622).
Function	domain:Alpha-helical parts of the C-terminal intracellular region mediate heterodimeric interaction with GABA-B receptor 1.,function:Receptor for GABA. The activity of this receptor is mediated by G-proteins that inhibit adenylyl cyclase activity, stimulates phospholipase A2, activates potassium channels, inactivates voltage-dependent calcium-channels and modulates inositol phospholipids hydrolysis. Plays a critical role in the fine-tuning of inhibitory synaptic transmission. Pre-synaptic GABA-B-R inhibit neurotransmitter release by down-regulating high-voltage activated calcium channels, whereas postsynaptic GABA-B-R decrease neuronal excitability by activating a prominent inwardly rectifying potassium (Kir) conductance that underlies the late inhibitory postsynaptic potentials. Not only implicated in synaptic inhibition but also in hippocampal long-term potentiation, slow wave sleep, m
Background	The multi-pass membrane protein encoded by this gene belongs to the G-protein coupled receptor 3 family and GABA-B receptor subfamily. The GABA-B receptors inhibit neuronal activity through G protein-coupled second-messenger systems, which regulate the release of neurotransmitters, and the activity of ion channels and adenylyl cyclase. This receptor subunit forms an active heterodimeric complex with GABA-B receptor subunit 1, neither of which is effective on its own. Allelic variants of this gene have been associated with nicotine dependence.[provided by RefSeq, Jan 2010],
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of 1 MCF-7 treated with LPS, 2 MCF7,using primary antibody at 1:1000 dilution. Secondary antibody(catalog#:RS23920) was diluted at 1:10000