



ENaCδ Polyclonal Antibody


Catalog No	YP-Ab-16412
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB;IF;ELISA
Gene Name	SCNN1D
Immunogen	The antiserum was produced against synthesized peptide derived from human SCNN1D. AA range:411-460
Specificity	ENaCδ Polyclonal Antibody detects endogenous levels of ENaCδ protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in other applications.
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Amiloride sensitive epithelial sodium channel gamma subunit;Amiloride sensitive sodium channel subunit gamma;Amiloride-sensitive sodium channel subunit gamma;BESC3;ENaC gamma subunit;ENaCG;ENaCgamma;Epithelial Na(+) channel subunit gamma;Epithelial Na+ channel subunit gamma;Gamma ENaC; Gamma NaCH;Gamma-ENaC;Gamma-NaCH;Nonvoltage gated sodium channel 1 subunit gamma;Nonvoltage-gated sodium channel 1 subunit gamma;PHA 1;PHA1;SCNEG;SCNN 1G;SCNN1G;SCNNG_HUMAN;Sodium channel non voltage gated 1 gamma subunit;Sodium channel nonvoltage gated 1 gamma;ENaCδ;SCNN1D;Amiloride-sensitive sodium channel subunit delta;DNACH;Delta-NaCH;Epithelial Na(+) channel subunit delta;Delta-ENaC;ENaCD; Nonvoltage-gated sodium channel 1 subunit delta;SCNED
Observed Band	70kD
Calculated Molecular Weight	70kD
Cell Pathway	Cell membrane ; Multi-pass membrane protein .
Tissue Specificity	Brain,Hippocampus,Lung,Nasal epithelium,Testis,




Function	function:Sodium permeable non-voltage-sensitive ion channel inhibited by the diuretic amiloride. Mediates the electrodiffusion of the luminal sodium (and water, which follows osmotically) through the apical membrane of epithelial cells. Controls the reabsorption of sodium in kidney, colon, lung and sweat glands. Also plays a role in taste perception.,similarity:Belongs to the amiloride-sensitive sodium channel family.,subunit:Heterotetramer of two alpha, one beta and one gamma subunit. A delta subunit can replace the alpha subunit.,
Background	function:Sodium permeable non-voltage-sensitive ion channel inhibited by the diuretic amiloride. Mediates the electrodiffusion of the luminal sodium (and water, which follows osmotically) through the apical membrane of epithelial cells. Controls the reabsorption of sodium in kidney, colon, lung and sweat glands. Also plays a role in taste perception.,similarity:Belongs to the amiloride-sensitive sodium channel family.,subunit:Heterotetramer of two alpha, one beta and one gamma subunit. A delta subunit can replace the alpha subunit.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.




Products Images



Immunofluorescence analysis of A549 cells, using SCNN1D Antibody. The picture on the right is blocked with the synthesized peptide.



Western blot analysis of lysates from RAW264.7 cells, using SCNN1D Antibody. The lane on the right is blocked with the synthesized peptide.



Western blot analysis of the lysates from HT-29 cells using SCNN1D antibody.