



MaxiK β 2 Polyclonal Antibody

Catalog No	YP-Ab-16465
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB;IHC;IF;ELISA
Gene Name	KCNMB2
Immunogen	The antiserum was produced against synthesized peptide derived from human KCNMB2. AA range:151-200
Specificity	MaxiK β 2 Polyclonal Antibody detects endogenous levels of MaxiK β 2 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB: 1/500 - 1/2000. IHC: 1/100 - 1/300. ELISA: 1/20000.. IF 1:50-200
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	2749B16Rik;31131N4Rik;BK channel subunit beta 2;BK channel subunit beta-2; BKbeta2;Calcium activated potassium channel subfamily M subunit beta 2; Calcium activated potassium channel subunit beta 2;Calcium-activated potassium channel;Calcium-activated potassium channel subunit beta-2; Charybdotoxin receptor subunit beta 2;Charybdotoxin receptor subunit beta-2; Hbeta2;Hbeta3;K(VCA)beta-2;K(VCA)beta2;KCMB2_HUMAN;KCNMB 2; Kcnmb2;Large conductance Ca ²⁺ activated K ⁺ channel beta2 subunit;Large conductance calcium activated potassium channel beta 2 subunit;Maxi K channel subunit beta 2;Maxi K channel subunit beta-2;MaxiK channel beta 2 subunit;MGC22431;MGC57945;Potassium large conductance calcium activated channel subfamily M beta member 2;Slo-beta-2;Slobeta2;subfamily M subunit beta-2;MaxiK β 2;Maxi K cha
Observed Band	30kD
Calculated Molecular Weight	30kD
Cell Pathway	Membrane; Multi-pass membrane protein.
Tissue Specificity	Expressed in kidney, heart and brain. Highly expressed in ovary. Expressed at low level in other tissues.



Function

domain:The ball and chain domain mediates the inactivation of KCNMA1. It occludes the conduction pathway of KCNMA1 channels, and comprises the pore-blocking ball domain (residues 1-17) and the chain domain (residues 20-45) linking it to the transmembrane segment. The ball domain is made up of a flexible N-terminus anchored at a well ordered loop-helix motif. The chain domain consists of a 4-turn helix with an unfolded linker at its C-terminus., function:Regulatory subunit of the calcium activated potassium KCNMA1 (maxiK) channel. Modulates the calcium sensitivity and gating kinetics of KCNMA1, thereby contributing to KCNMA1 channel diversity. Acts as a negative regulator that confers rapid and complete inactivation of KCNMA1 channel complex. May participate in KCNMA1 inactivation in chromaffin cells of the adrenal gland or in hippocampal CA1 neurons.,PTM:N-glycosylated.,similarity:Belongs

Background

MaxiK channels are large conductance, voltage and calcium-sensitive potassium channels which are fundamental to the control of smooth muscle tone and neuronal excitability. MaxiK channels can be formed by 2 subunits: the pore-forming alpha subunit and the modulatory beta subunit. The protein encoded by this gene is an auxiliary beta subunit which decreases the activation time of MaxiK alpha subunit currents. Alternative splicing results in multiple transcript variants of this gene. Additional variants are discussed in the literature, but their full length nature has not been described. [provided by RefSeq, Jul 2013],

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions


This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using MaxiKβ2 Polyclonal Antibody diluted at 1:500



Western blot analysis of lysates from Jurkat, COLO, and HepG2 cells, using KCNMB2 Antibody. The lane on the right is blocked with the synthesized peptide.



Western blot analysis of the lysates from HepG2 cells using KCNMB2 antibody.

