




Cdc50B Polyclonal Antibody


Catalog No	YP-Ab-16687
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB;ELISA
Gene Name	TMEM30B
Immunogen	The antiserum was produced against synthesized peptide derived from human TMEM30B. AA range:91-140
Specificity	Cdc50B Polyclonal Antibody detects endogenous levels of Cdc50B protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	91311B11Rik;CC5B_HUMAN;CDC5;S. cerevisiae;homolog of;B;CDC5B;Cell cycle control protein 5B;MGC126775;TMEM3B;Transmembrane protein 3B;Cdc50B;TMEM30B;Cell cycle control protein 50B;Transmembrane protein 30B
Observed Band	38kD
Calculated Molecular Weight	38kD
Cell Pathway	Cell membrane ; Multi-pass membrane protein .
Tissue Specificity	Heart,Lung,
Function	similarity:Belongs to the CDC50/LEM3 family.,
Background	similarity:Belongs to the CDC50/LEM3 family.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.




Products Images



Western Blot analysis of K562 cells using Cdc50B Polyclonal Antibody



Western blot analysis of lysates from K562 cells, using TMEM30B Antibody. The lane on the right is blocked with the synthesized peptide.



Western blot analysis of the lysates from HeLa cells using TMEM30B antibody.