




CNTD2 Polyclonal Antibody

Catalog No	YP-Ab-16715
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	IHC;IF;ELISA
Gene Name	CNTD2
Immunogen	The antiserum was produced against synthesized peptide derived from human CNTD2. AA range:51-100
Specificity	CNTD2 Polyclonal Antibody detects endogenous levels of CNTD2 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	IHC: 1/100 - 1/300. ELISA: 1/40000.. IF 1:50-200
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Homo sapiens (Human);Cyclin-P;Cyclin P;CNTD2;Cyclin N-terminal domain-containing protein 2
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Nucleus .
Tissue Specificity	Expression is induced in lung and colon cancer cells.
Function	similarity:Belongs to the cyclin family.,similarity:Contains 1 cyclin N-terminal domain.,
Background	similarity:Belongs to the cyclin family.,similarity:Contains 1 cyclin N-terminal domain.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Immunohistochemistry analysis of paraffin-embedded human breast carcinoma, using CNTD2 Antibody. The picture on the right is blocked with the synthesized peptide.