



KIAA0100 Polyclonal Antibody

Catalog No	YP-Ab-16743
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	IHC;IF;ELISA
Gene Name	KIAA0100
Immunogen	The antiserum was produced against synthesized peptide derived from human KIAA0100. AA range:681-730
Specificity	KIAA0100 Polyclonal Antibody detects endogenous levels of KIAA0100 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in other applications.
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Antigen MLAA-22;BCOX;BCOX1;Breast cancer-overexpressed gene 1 protein; CT11;DKFZp686M843;K1_HUMAN;KIAA1;MGC111488;MGC134981;UPF378 protein KIAA1;KIAA0100;UPF0378 protein KIAA0100
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Secreted .
Tissue Specificity	Expressed in pancreas, placenta and up-regulated in breast carcinoma epithelial cells, ductal in situ carcinoma (DCIS), invasive breast carcinoma (IBC) and metastatic breast carcinoma cells (MET).
Function	sequence caution:Contaminating sequence.,sequence caution:Sequencing errors.,similarity:Belongs to the UPF0378 family.,tissue specificity:Expressed in pancreas, placenta and up-regulated in breast carcinoma epithelial cells, ductal in situ carcinoma (DCIS), invasive breast carcinoma (IBC) and metastatic breast carcinoma cells (MET).,



Background

sequence caution:Contaminating sequence.,sequence caution:Sequencing errors.,similarity:Belongs to the UPF0378 family.,tissue specificity:Expressed in pancreas, placenta and up-regulated in breast carcinoma epithelial cells, ductal in situ carcinoma (DCIS), invasive breast carcinoma (IBC) and metastatic breast carcinoma cells (MET).,


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Immunohistochemistry analysis of paraffin-embedded human breast carcinoma, using K0100 Antibody. The picture on the right is blocked with the synthesized peptide.