



p18 Polyclonal Antibody

Catalog No	YP-Ab-16760
Isotype	IgG
Reactivity	Human;Mouse;Monkey
Applications	WB;IHC;IF;ELISA
Gene Name	CDKN2C
Immunogen	The antiserum was produced against synthesized peptide derived from human p18 INK. AA range:111-160
Specificity	p18 Polyclonal Antibody detects endogenous levels of p18 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB: 1/500 - 1/2000. IHC: 1/100 - 1/300. ELISA: 1/20000.. IF 1:50-200
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	UBE2I;C358B7.1;E2;EC:2.3.2.-;P18;CDKN2C;Cyclin-dependent kinase 4 inhibitor C;CDKN6;Cyclin-dependent kinase 6 inhibitor;p18-INK4c;p18-INK6
Observed Band	18kD
Calculated Molecular Weight	18kD
Cell Pathway	nucleus,cytoplasm,cytosol,
Tissue Specificity	Highest levels found in skeletal muscle. Also found in pancreas and heart.
Function	disease:Defects in CDKN2C are involved in tumor formation.,function:Interacts strongly with CDK6, weakly with CDK4. Inhibits cell growth and proliferation with a correlated dependence on endogenous retinoblastoma protein RB., similarity:Belongs to the CDKN2 cyclin-dependent kinase inhibitor family., similarity:Contains 4 ANK repeats.,subunit:Heterodimer of p18 with CDK6., tissue specificity:Highest levels found in skeletal muscle. Also found in pancreas and heart.,
Background	The protein encoded by this gene is a member of the INK4 family of cyclin-dependent kinase inhibitors. This protein has been shown to interact with CDK4 or CDK6, and prevent the activation of the CDK kinases, thus function as a cell



growth regulator that controls cell cycle G1 progression. Ectopic expression of this gene was shown to suppress the growth of human cells in a manner that appears to correlate with the presence of a wild-type RB1 function. Studies in the knockout mice suggested the roles of this gene in regulating spermatogenesis, as well as in suppressing tumorigenesis. Two alternatively spliced transcript variants of this gene, which encode an identical protein, have been reported. [provided by RefSeq, Jul 2008],


matters needing attention

Avoid repeated freezing and thawing!


Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.


Products Images



Western Blot analysis of various cells using p18 Polyclonal Antibody



Immunohistochemistry analysis of paraffin-embedded human breast carcinoma tissue, using p18 INK Antibody. The picture on the right is blocked with the synthesized peptide.



Western blot analysis of lysates from COS7 cells, using p18 INK Antibody. The lane on the right is blocked with the synthesized peptide.