



Actinin- α 1/2/3/4 Polyclonal Antibody

Catalog No	YP-Ab-16923
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB;IHC;IF;ELISA
Gene Name	ACTN1/2/3/4
Immunogen	The antiserum was produced against synthesized peptide derived from human ACTN1/2/3/4. AA range:21-70
Specificity	Actinin- α 1/2/3/4 Polyclonal Antibody detects endogenous levels of Actinin- α 1/2/3/4 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB: 1/500 - 1/2000. IHC: 1/100 - 1/300. ELISA: 1/20000.. IF 1:50-200
Concentration	1 mg/ml
Purity	$\geq 90\%$
Storage Stability	-20°C/1 year
Synonyms	actinin 1 smooth muscle;Actinin alpha 1;actinin;alpha 1;ACTN 1;Actn1; ACTN1_HUMAN;Alpha Actinin 1;Alpha actinin cytoskeletal isoform;Alpha-actinin cytoskeletal isoform;Alpha-actinin-1;BDPLT15;F actin cross linking protein;F-actin cross-linking protein;FLJ4884;FLJ54432;Non muscle alpha actinin 1;Non-muscle alpha-actinin-1;Actin binding protein;Actinin alpha 2;ACTN 2;ACTN2; ACTN2_HUMAN;Alpha actinin 2;Alpha actinin skeletal muscle;Alpha actinin skeletal muscle isoform 2;Alpha-actinin skeletal muscle isoform 2;Alpha-actinin-2;CMD1AA;Actinin alpha 3;alpha;ACTN 3;ACTN3;ACTN3_HUMAN;Alpha actinin skeletal muscle isoform 3;Alpha-Actinin 3;Alpha-actinin skeletal muscle isoform 3;Alpha-actinin-3;MGC1172;MGC1175;OTTHUMP21949; OTTHUMP21954;actinin 4;Actinin alpha 4;actinin4;ACTN 4;ACTN4; ACTN4_HUMAN;Alpha-actinin-4;DKFZp686K23158;Focal segmental glomerulosclerosis 1;FSGS 1;FSGS;FSGS1;Non muscle alpha actinin 4;Non-muscle alpha-actinin 4;Actinin- α 1/2/3/4;ACTN1/2/3/4;Alpha-actinin-1/2/3/4;A
Observed Band	75kD
Calculated Molecular Weight	75kD
Cell Pathway	Cytoplasm, cytoskeleton . Cytoplasm, myofibril, sarcomere, Z line . Cell membrane . Cell junction . Cell projection, ruffle . Colocalizes with MYOZ2 and PPP3CA at the Z-line of heart and skeletal muscle. Colocalizes with PSD in




membrane ruffles and central reticular structures. .

Tissue Specificity	Colon,Epithelium,Hepatocyte,Neuroblastoma,Placenta,Platelet,Skin,Tumor of tongue,
Function	function:F-actin cross-linking protein which is thought to anchor actin to a variety of intracellular structures. This is a bundling protein.,similarity:Belongs to the alpha-actinin family.,similarity:Contains 1 actin-binding domain.,similarity:Contains 2 CH (calponin-homology) domains.,similarity:Contains 2 EF-hand domains.,similarity:Contains 4 spectrin repeats.,subcellular location:Colocalizes with MYOZ2 and PPP3CA at the Z-line of heart and skeletal muscle.,subunit:Homodimer; antiparallel. Interacts with DDN, MYOZ2, PDLIM2, TTID and LPP.,
Background	Alpha actinins belong to the spectrin gene superfamily which represents a diverse group of cytoskeletal proteins, including the alpha and beta spectrins and dystrophins. Alpha actinin is an actin-binding protein with multiple roles in different cell types. In nonmuscle cells, the cytoskeletal isoform is found along microfilament bundles and adherens-type junctions, where it is involved in binding actin to the membrane. In contrast, skeletal, cardiac, and smooth muscle isoforms are localized to the Z-disc and analogous dense bodies, where they help anchor the myofibrillar actin filaments. This gene encodes a nonmuscle, cytoskeletal, alpha actinin isoform and maps to the same site as the structurally similar erythroid beta spectrin gene. Three transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2008],
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.




Products Images




Western Blot analysis of various cells using Actinin- α 1/2/3/4 Polyclonal Antibody diluted at 1:1000



Western Blot analysis of A549 cells using Actinin- α 1/2/3/4 Polyclonal Antibody diluted at 1:1000



Immunohistochemical analysis of paraffin-embedded Human brain. Antibody was diluted at 1:100(4° overnight). High-pressure and temperature Tris-EDTA, pH8.0 was used for antigen retrieval. Negative control (right) obtained from antibody was pre-absorbed by immunogen peptide.



Western blot analysis of lysates from 293, HeLa, and A549 cells, using ACTN1/2/3/4 Antibody. The lane on the right is blocked with the synthesized peptide.