



CEACAM21 Polyclonal Antibody


Catalog No	YP-Ab-16949
Isotype	IgG
Reactivity	Human
Applications	WB;ELISA
Gene Name	CEACAM21
Immunogen	The antiserum was produced against synthesized peptide derived from human CEACAM21. AA range:171-220
Specificity	CEACAM21 Polyclonal Antibody detects endogenous levels of CEACAM21 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Carcinoembryonic antigen-related cell adhesion molecule 21;UNQ3098/PRO10075;CEACAM21
Observed Band	32kD
Calculated Molecular Weight	32kD
Cell Pathway	Membrane ; Single-pass type I membrane protein .
Tissue Specificity	Testis,
Function	similarity:Belongs to the immunoglobulin superfamily. CEA family., similarity:Contains 1 Ig-like C2-type (immunoglobulin-like) domain.,
Background	similarity:Belongs to the immunoglobulin superfamily. CEA family., similarity:Contains 1 Ig-like C2-type (immunoglobulin-like) domain.,
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of the lysates from HT-29 cells using CEACAM21 antibody.