



COL19A1 Polyclonal Antibody

Catalog No	YP-Ab-16976
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB;IHC;IF;ELISA
Gene Name	COL19A1
Immunogen	The antiserum was produced against synthesized peptide derived from human Collagen XIX alpha1. AA range:421-470
Specificity	COL19A1 Polyclonal Antibody detects endogenous levels of COL19A1 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	COJA1_HUMAN;COL19A1;COL9A1L;Collagen alpha-1(XIX) chain;Collagen alpha-1(Y) chain;D6S228E
Observed Band	120kD
Calculated Molecular Weight	120kD
Cell Pathway	Secreted, extracellular space, extracellular matrix .
Tissue Specificity	Localized to vascular, neuronal, mesenchymal, and some epithelial basement membrane zones in umbilical cord.
Function	domain:The numerous interruptions in the triple helix may make this molecule either elastic or flexible.,function:May act as a cross-bridge between fibrils and other extracellular matrix molecules. Involved in skeletal myogenesis in the developing esophagus. May play a role in organization of the pericellular matrix or the sphinteric smooth muscle.,PTM:Prolines at the third position of the tripeptide repeating unit (G-X-Y) are hydroxylated in some or all of the chains.,similarity:Belongs to the fibril-associated collagens with interrupted helices (FACIT) family.,similarity:Contains 1 TSP N-terminal (TSPN) domain.,subunit:Oligomer; disulfide-linked.,tissue specificity:Localized to vascular,



neuronal, mesenchymal, and some epithelial basement membrane zones in umbilical cord.,

Background

This gene encodes the alpha chain of type XIX collagen, a member of the FACIT collagen family (fibril-associated collagens with interrupted helices). Although the function of this collagen is not known, other members of this collagen family are found in association with fibril-forming collagens such as type I and II, and serve to maintain the integrity of the extracellular matrix. The transcript produced from this gene has an unusually large 3' UTR which has not been completely sequenced. [provided by RefSeq, Jul 2008],


matters needing attention

Avoid repeated freezing and thawing!


Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Immunofluorescence analysis of A549 cells, using Collagen XIX alpha1 Antibody. The picture on the right is blocked with the synthesized peptide.



Immunohistochemistry analysis of paraffin-embedded human thymus gland tissue, using Collagen XIX alpha1 Antibody. The picture on the right is blocked with the synthesized peptide.