



PCYT2 rabbit pAb

Catalog No	YP-Ab-17123
Isotype	IgG
Reactivity	Human, Mouse, Rat
Applications	WB 1:500 - 1:2000 Polyclonal Antibody
Immunogen	Recombinant fusion protein of human PCYT2 (NP_002852.1).
matters needing attention	Avoid repeated freezing and thawing!
Formulation	
Host	Rabbit
Purification	Affinity purification
Dilution	
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	type(□□): Polyclonal Antibody
Observed Band	44KD
Calculated Molecular Weight	
Cell Pathway	
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.
Disclaimer note	The observed molecular weight of the protein may vary from the listed predicted molecular weight due to post translational modifications, post translation cleavages, relative charges, and other experimental factors.
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
