



TAP-Tag Monoclonal Antibody(Mix), HRP Conjugated

Catalog No	YP-Ab-17205
Isotype	IgG
Reactivity	Species independent
Applications	WB
Gene Name	
Immunogen	
Specificity	TAP Tag Monoclonal Antibody(Mix) HRP Conjugated, specially designed for your Western blot analysis.
Formulation	Liquid in PBS, pH 7.4, containing 0.02% sodium azide as preservative and 50% Glycerol.
Source	Monoclonal, Mouse IgG
Purification	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Dilution	Optimal working dilutions should be determined experimentally by the investigator. Suggested starting dilutions are as follows: WB 1:5000-10000.
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	TAP;TAP tag;TAP-tag
Observed Band	75kDa (Recombinant TAP protein)
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

