



TurboGFP Mouse mAb(Mix-mA)

Catalog No	YP-Ab-17214
Isotype	IgG
Reactivity	Species independent
Applications	WB
Gene Name	TurboGFP
Immunogen	Synthesized peptide derived from human TurboGFP
Specificity	This antibody detects endogenous levels of TurboGFP at Species independent
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Mouse, monoclonal
Purification	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	TurboGFP
Observed Band	26kD
Calculated Molecular Weight	26kD
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
