



Zip9 Mouse mAb(Mix-mA)

Catalog No	YP-Ab-17216
Isotype	IgG
Reactivity	Human, Mouse,Rat
Applications	WB,IHC
Gene Name	SLC39A9 ZIP9 UNQ714/PRO1377
Immunogen	Synthesized peptide derived from human Zip9
Specificity	This antibody detects endogenous levels of Zip9 at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Mouse,monoclonal
Purification	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 IHC 1:50-200
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	SLC39A9;Zinc transporter ZIP9;ZIP 9;ZIP9;Zrt and Irt like protein 9;SLC39A9 ZIP9 UNQ714/PRO1377;Zinc transporter ZIP9 (Solute carrier family 39 member 9) (Zrt- and Irt-like protein 9) (ZIP-9)
Observed Band	25 kDa
Calculated Molecular Weight	307 aa, 32 kDa
Cell Pathway	Membrane ; Multi-pass membrane protein .
Tissue Specificity	
Function	function:May act as a zinc-influx transporter.,similarity:Belongs to the ZIP transporter (TC 2.A.5) family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
