




Mitoregulin rabbit pAb

Catalog No	YP-Ab-17229
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	LINC00116 NCRNA00116 MTLN
Immunogen	Synthesized peptide derived from human YB039 AA:20-40
Specificity	This antibody detects endogenous levels of YB039 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	YB039;LINC00116 NCRNA00116 MTLN;YB039 Mitoregulin
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Mitochondrion inner membrane ; Single-pass membrane protein . Preferentially binds to cardiolipin relative to other common cell membrane lipids. .
Tissue Specificity	
Function	caution:It is uncertain whether Met-1 or Met-83 is the initiator.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from MCF-7 cells,
primary antibody was diluted at 1:1000, 4° over night