



UBL4A rabbit pAb

Catalog No	YP-Ab-17238
Isotype	IgG
Reactivity	Human, Mouse,Rat
Applications	IHC,WB
Gene Name	UBL4A DXS254E GDX UBL4
Immunogen	Synthesized peptide derived from human N-teral UBL4A
Specificity	This antibody detects endogenous levels of UBL4A at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Rabbit,polyclonal
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000 IHC 1:50-200
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Ubiquitin-like protein 4A;ubiquitin like 4A;GDX;G6PD;DXS254E;UBL4A;UBL4A DXS254E GDX UBL4;Ubiquitin-like protein 4A (Ubiquitin-like protein GDX)
Observed Band	17 kDa
Calculated Molecular Weight	18 kDa
Cell Pathway	Cytoplasm, cytosol . Nucleus .
Tissue Specificity	
Function	As part of a cytosolic protein quality control complex, the BAG6/BAT3 complex, maintains misfolded and hydrophobic patches-containing proteins in a soluble state and participates in their proper delivery to the endoplasmic reticulum or alternatively can promote their sorting to the proteasome where they undergo degradation . The BAG6/BAT3 complex is involved in the post-translational delivery of tail-anchored/type II transmembrane proteins to the endoplasmic reticulum membrane. Recruited to ribosomes, it interacts with the transmembrane region of newly synthesized tail-anchored proteins and together with SGTA and ASNA1 mediates their delivery to the endoplasmic reticulum . Client proteins that cannot be properly delivered to the endoplasmic reticulum are ubiquitinated and sorted to the proteasome . Similarly, the BAG6/BAT3 complex also functions as a sorting platform for proteins of the

Background

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images