



LOX12 rabbit pAb

Catalog No	YP-Ab-17240
Isotype	IgG
Reactivity	Human, Mouse
Applications	IHC, WB
Gene Name	ALOX12 LOG12
Immunogen	Synthesized peptide derived from human N-terminal LOX12
Specificity	This antibody detects endogenous levels of LOX12 at Human, Mouse
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Rabbit, polyclonal
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000 IHC 1:50-200
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	LOG12; 12-LOX; 12S-LOX; ALOX12; ALOX12 LOG12; Arachidonate 12-lipoxygenase; 12S-type (12S-LOX) (12S-lipoxygenase) (EC 1.13.11.31) (Platelet-type lipoxygenase 12)
Observed Band	76 kDa
Calculated Molecular Weight	73kD
Cell Pathway	Cytoplasm, cytosol. Membrane. Membrane association is stimulated by EGF.
Tissue Specificity	Expressed in vascular smooth muscle cells.
Function	Catalyzes the regio and stereo-specific incorporation of molecular oxygen into free and esterified polyunsaturated fatty acids generating lipid hydroperoxides that can be further reduced to the corresponding hydroxy species. Mainly converts arachidonate ((5Z,8Z,11Z,14Z)-eicosatetraenoate) to the specific bioactive lipid (12S)-hydroperoxyeicosatetraenoate/(12S)-HPETE. Through the production of bioactive lipids like (12S)-HPETE it regulates different biological processes including platelet activation. It can also catalyze the epoxidation of double bonds of polyunsaturated fatty acids such as (14S)-hydroperoxy-docosahexaenoate/(14S)-HPDHA resulting in the formation of (13S,14S)-epoxy-DHA. Furthermore, it may participate in the sequential oxidations of DHA ((4Z,7Z,10Z,13Z,16Z,19Z)-docosahexaenoate) to generate specialized pro-resolving

mediators (SPMs) like resolvin D5 ((7S,17S)-diHPDHA

Background

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images