



TFAP4 rabbit pAb

Catalog No	YP-Ab-17249
Isotype	IgG
Reactivity	Human
Applications	IHC, WB
Gene Name	TFAP4 BHLHC41
Immunogen	Synthesized peptide derived from human N-terminal TFAP4
Specificity	This antibody detects endogenous levels of TFAP4 at Human
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Rabbit, polyclonal
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000 IHC 1:50-200
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	AP 4; bHLHc41; TFAP4; Transcription factor AP 4; TFAP4 BHLHC41; Transcription factor AP-4 (Activating enhancer-binding protein 4) (Class C basic helix-loop-helix protein 41) (bHLHc41)
Observed Band	39-42 kDa
Calculated Molecular Weight	338 aa, 39 kDa
Cell Pathway	Nucleus.
Tissue Specificity	
Function	Transcription factor that activates both viral and cellular genes by binding to the symmetrical DNA sequence 5'-CAGCTG-3'.
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
