



PTPRR rabbit pAb

Catalog No	YP-Ab-17251
Isotype	IgG
Reactivity	Human, Mouse,Rat
Applications	IHC,WB
Gene Name	PTPRR ECPTP PTPRQ
Immunogen	Synthesized peptide derived from human N-ternal PTPRR
Specificity	This antibody detects endogenous levels of PTPRR at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Rabbit,polyclonal
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000 IHC 1:50-200
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Ch-1PTPase;EC PTP;EC:3.1.3.48;ECPTP;NC-PTPCOM1;PTPRR;PTPRR ECPTP PTPRQ;Receptor-type tyrosine-protein phosphatase R (R-PTP-R) (EC 3.1.3.48) (Ch-1PTPase) (NC-PTPCOM1) (Protein-tyrosine phosphatase PCPTP1)
Observed Band	32 kDa, 47 kDa, 65 kDa
Calculated Molecular Weight	657 aa, 74 kDa
Cell Pathway	[Isoform Alpha]: Cell membrane; Single-pass type I membrane protein.; [Isoform Delta]: Cytoplasm, perinuclear region. Locates to the perinuclear areas within the cytoplasm.; [Isoform Gamma]: Cytoplasm, perinuclear region. Locates to the perinuclear areas within the cytoplasm.
Tissue Specificity	Expressed in brain, placenta, small intestine, stomach, uterus and weakly in the prostate. Isoform alpha has been observed only in the brain. Isoform gamma is expressed in brain, placenta and uterus. Isoform delta is expressed in brain, kidney, placenta, prostate, small intestine and uterus.
Function	Sequesters mitogen-activated protein kinases (MAPKs) such as MAPK1, MAPK3 and MAPK14 in the cytoplasm in an inactive form. The MAPKs bind to a dephosphorylated kinase interacting motif, phosphorylation of which by the protein kinase A complex releases the MAPKs for activation and translocation

into the nucleus (By similarity).

Background

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images