



FBX16 rabbit pAb

Catalog No	YP-Ab-17262
Isotype	IgG
Reactivity	Human, Mouse
Applications	IHC, WB
Gene Name	FBXO16 FBX16
Immunogen	Synthesized peptide derived from human C-terminal FBX16
Specificity	This antibody detects endogenous levels of FBX16 at Human, Mouse
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Rabbit, polyclonal
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000 IHC 1:50-200
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	F box only protein 16; F box protein 16; FBX16; FBXO16; FBXO16 FBX16; F-box only protein 16
Observed Band	35 kDa
Calculated Molecular Weight	292 aa, 35 kDa
Cell Pathway	
Tissue Specificity	Expressed in heart, spleen and colon.
Function	Probably recognizes and binds to some phosphorylated proteins and promotes their ubiquitination and degradation.
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
