



S43A3 rabbit pAb

Catalog No	YP-Ab-17264
Isotype	IgG
Reactivity	Human, Mouse
Applications	IHC, WB
Gene Name	SLC43A3 PSEC0252
Immunogen	Synthesized peptide derived from human C-terminal S43A3
Specificity	This antibody detects endogenous levels of S43A3 at Human, Mouse
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Rabbit, polyclonal
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000 IHC 1:50-200
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	SLC43A3; SLC43A3 PSEC0252; Solute carrier family 43 member 3 (Protein FOAP-13)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Basolateral cell membrane ; Multi-pass membrane protein .
Tissue Specificity	Widely expressed with highest levels in the liver and lung, followed by the pancreas (PubMed:26455426). Highly expressed in macrophages (Ref.1).
Function	Sodium-independent purine-selective nucleobase transporter which mediates the equilibrative transport of extracellular purine nucleobases such as adenine, guanine and hypoxanthine . May regulate fatty acid (FA) transport in adipocytes, acting as a positive regulator of FA efflux and as a negative regulator of FA uptake (By similarity). ; [Isoform 1]: Sodium-independent purine-selective nucleobase transporter which mediates the equilibrative transport of extracellular purine nucleobase adenine . Mediates the influx and efflux of the purine nucleobase analog drug 6-mercaptopurine across the membrane . ; [Isoform 2]: Sodium-independent purine-selective nucleobase transporter which mediates the equilibrative transport of extracellular purine nucleobase adenine . Mediates the influx and efflux of the purine nucleobase analog drug 6-



mercaptopurine across the membrane .

Background

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images