



GCSH rabbit pAb

Catalog No	YP-Ab-17289
Isotype	IgG
Reactivity	Human, Mouse,Rat
Applications	IHC,WB
Gene Name	GCSH
Immunogen	Synthesized peptide derived from human C-ternal GCSH
Specificity	This antibody detects endogenous levels of GCSH at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Rabbit,polyclonal
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000 IHC 1:50-200
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	NKH;MMDS7;Glycine cleavage system H protein;mitochondrial;GCE;GCSH
Observed Band	15 kDa
Calculated Molecular Weight	19 kDa
Cell Pathway	Mitochondrion .
Tissue Specificity	
Function	The glycine cleavage system catalyzes the degradation of glycine. The H protein (GCSH) shuttles the methylamine group of glycine from the P protein (GLDC) to the T protein (GCST).
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

