



Anti-Covid-19 Spike mAb IgG □ S04FC)

Catalog No	YP-Ab-17292
Isotype	IgG
Reactivity	Covid-19
Applications	ELISA ECL Immunogold
Gene Name	SARS-CoV-2 S, 2019-nCoV S, Covid-19 S
Immunogen	The recombinant SARS-CoV-2 (2019-nCoV) Spike Protein S
Specificity	This antibody detects SARS-CoV-2 (2019-nCoV) Spike Protein S1 doman
Formulation	PBS
Source	Humanized mAb IgG
Purification	The antibody was affinity-purified from 293T expression by affinity-chromatography using specific immunogen.
Dilution	ELISA 1:5000-1:100000, WB 1:1000-1:10000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Anti-Covid-19 Spike mAb IgG(S04FC);SARS-CoV-2 S;2019-nCoV S;Covid-19 S; SARS-CoV-2 Spike;2019-nCoV Spike;Covid-19 Spike
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

