



Anti-Covid-19 N protein mAb IgM□ S06FC)

Catalog No	YP-Ab-17296
Isotype	IgG
Reactivity	Covid-19
Applications	ELISA ECL Immunogold
Gene Name	SARS-CoV-2 N, 2019-nCoV N, Covid-19 N
Immunogen	The recombinant SARS-CoV-2 (2019-nCoV) N Protein
Specificity	This antibody detects SARS-CoV-2 (2019-nCoV) Nucleoprotein
Formulation	PBS
Source	Humanized mAb IgM
Purification	The antibody was affinity-purified from 293T expression by affinity-chromatography using specific immunogen.
Dilution	ELISA 1:5000-1:100000, WB 1:1000-1:10000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Anti-Covid-19 N protein mAb IgM(S06FC);SARS-CoV-2 N;2019-nCoV N;Covid-19 N;SARS-CoV-2 Nucleoprotein;2019-nCoV Nucleoprotein;Covid-19 Nucleoprotein
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
