



Glypican 3 Rabbit mAb

Catalog No	YP-Ab-17684
Isotype	IgG
Reactivity	Human
Applications	WB,FC
Gene Name	GPC3
Research Field	Stem Cells
Product Categories	Primary antibody
Host	Rabbit
Clonality	Monoclonal Antibody
Clonality No.	R02-8C8
Immunogen	A synthesized peptide derived from human Glypican 3 Purification Affinity Chromatography Conjugation Unconjugated Modification Unmodified Form Liquid
Buffer System	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Dilution	WB: 1/500-1/1000 FC: 1/50-1/100
Storage Stability	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Synonyms	GPC3;DGSX;Glypican-3 alpha subunit;Glypican-3 beta subunit;GTR2 2;Dnmt3b; DNA methyltransferase HsaIIIB;OCI5;Glypican-3;GTR2-2;Intestinal protein OCI-5;MXR7
SwissProt ID	P51654
Gene ID	2719
Molecular Weight	580 aa, 66 kDa
Background	Cell surface proteoglycan that bears heparan sulfate. Inhibits the dipeptidyl peptidase activity of DPP4. May be involved in the suppression/modulation of growth in the predominantly mesodermal tissues and organs.
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

