



MTTP/MTP rabbit pAb

Catalog No	YP-Ab-17686
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Target ☐	MTTP
Gene Name	MTTP MTP
Immunogen	Synthesized peptide derived from human MTTP/MTP
Specificity ☐	This antibody detects endogenous levels of MTTP/MTP at Human, Mouse
Formulation ☐	Liquid in PBS containing 50% glycerol, and 0.3090% sodium azide.
Source ☐	Rabbit,polyclonal
Expression ☐	Liver and small intestine. Also found in ovary, testis and kidney.
Purification ☐	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration ☐	1 mg/ml
Storage Stability	-15°C to -25°C/1 year(Do not lower than -25°C)
Molecular Weight(Da) ☐	98kD
Function ☐	Catalyzes the transport of triglyceride, cholesteryl ester, and phospholipid between phospholipid surfaces . Required for the assembly and secretion of plasma lipoproteins that contain apolipoprotein B . May be involved in regulating cholesteryl ester biosynthesis in cells that produce lipoproteins (By similarity).
Subcellular Location ☐	Endoplasmic reticulum . Golgi apparatus . Colocalizes with P4HB/PDI in the endoplasmic reticulum (PubMed:23475612, PubMed:26224785).
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images