



CCR2 Rabbit pAb

Catalog No	YP-Ab-17695
Isotype	IgG
Reactivity	Human, Mouse
Applications	WB, IHC-P
Gene Name	CCR2
Research Field	Immunology
Host	Rabbit
Immunogen	Synthetic peptide of human CCR2
Buffer System	pH 7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification	Affinity Purified Conjugation Unconjugated
Dilution	WB: 1/500-1/1000 IHC: 1/50-1/100
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Synonyms	CKR2; CCR-2; CCR2A; CCR2B; CD192; CKR2A; CKR2B; CMKBR2; MCP-1-R; CC-CKR-2; CCR2; Cytokeratin (Pan); C-C chemokine receptor type 2; C-C CKR-2; Monocyte chemoattractant protein 1 receptor; CD antigen CD192
Molecular Weight	Calculated MW: 42 kDa; Observed MW: 42 kDa
Background	Receptor for the MCP-1, MCP-3 and MCP-4 chemokines. Transduces a signal by increasing the intracellular calcium ions level. Alternative coreceptor with CD4 for HIV-1 infection.
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

