



KCTD12 rabbitt pAb

□□	YP-Ab-17696 IgG Human;Mouse
□□	WB
□□□	KCTD12
□□□□□	KCTD12 C13orf2 KIAA1778 PFET1
□□□□	Synthesized peptide derived from human KCTD12
□□□□	This antibody detects endogenous levels of KCTD12 at Human, Mouse Liquid in PBS containing 50% glycerol, and 0.305% sodium azide.
□□□	
□□□	Polyclonal, Rabbit,IgG
□□□	WB 1:500-2000
□□□□□	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
□□	1 mg/ml
□□	≥90%
□□□	-15°C to -25°C/1 year(Do not lower than -25°C)
□□□□	36 kDa
□□□	Auxiliary subunit of GABA-B receptors that determine the pharmacology and kinetics of the receptor response. Increases agonist potency and markedly alter the G-protein signaling of the receptors by accelerating onset and promoting desensitization (By similarity).
□□□□□	Cell junction, synapse, presynaptic cell membrane. Cell junction, synapse, postsynaptic cell membrane .
□□□□□	Present in a variety of fetal organs, with highest expression levels in the cochlea and brain and, in stark contrast, is detected only at extremely low levels in adult organs, such as brain and lung.
□□□□□□	Avoid repeated freezing and thawing!
□□□□	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

