



His Tag Mouse mAb (HRP)

Catalog No	YP-Ab-17699
Isotype	IgG
Alternative Names	6 His epitope tag; Hexa His tag; HHHHHH epitope tag; HHHHHH tag; His tag
Research Field	Tags & Cell Markers
Product Categories	Tag antibody
Host	Mouse
Reactivity	Species-independent
Application	WB
Dilution Ratio	WB: 1/1000-1/10000
Molecular Weight	Refer to figures
Clonality	Monoclonal Antibody
Immunogen	This His Tag antibody is generated from a mouse immunized with recombinant protein.
Purification	Affinity Purified
Conjugation	HRP
Modification	Unmodified
Form	Liquid
Buffer System	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Background	A polyhistidine-tag is an amino acid motif in proteins that consists of at least five histidine (His) residues, often at the N- or C-terminus of the protein. It is also known as hexa histidine-tag, 6xHis-tag, and by the trademarked name His-tag. Polyhistidine-tags are often used for affinity purification of polyhistidine-tagged recombinant proteins expressed in Escherichia coli and other prokaryotic expression systems.
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

