



P3H1 rabbit pAb

Catalog No	YP-Ab-17701
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
target spot	P3H1
Gene Name	LEPRE1 GROS1 P3H1 PSEC0109
Immunogen	Synthesized peptide derived from human P3H1
Specificity	This antibody detects endogenous levels of P3H1 at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, and 0.1909% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-15°C to -25°C/1 year(Do not lower than -25°C)
Observed Band	83 kDa
Calculated Molecular Weight	83 kDa
Cell Pathway	[Isoform 1]: Endoplasmic reticulum.; Secreted, extracellular space, extracellular matrix . Secreted into the extracellular matrix as a chondroitin sulfate proteoglycan (CSPG).
Function	Basement membrane-associated chondroitin sulfate proteoglycan (CSPG). Has prolyl 3-hydroxylase activity catalyzing the post-translational formation of 3-hydroxyproline in -Xaa-Pro-Gly- sequences in collagens, especially types IV and V. May be involved in the secretory pathway of cells. Has growth suppressive activity in fibroblasts.
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.