



CD168 Mouse mAb

Catalog No	YP-Ab-17714
Isotype	Mouse IgG1
Gene Name	CD168
Alternative Names	HMMR; IHABP; RHAMM
Research Field	Stem Cells
Product Categories	Primary antibody
Host	Mouse
Reactivity	Human
Application	FC,ELISA
Dilution Ratio	FC: 1/50-1/100 ELISA: 1/10000
Clonality	Monoclonal Antibody
Clonality No.	7J6-X6-M5
Immunogen	Purified recombinant fragment of human CD168 expressed in E. Coli.
Purification	Affinity Purified
Conjugation	Unconjugated
Modification	Unmodified
Form	Liquid
Buffer System	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Calculated Molecular Weight	
Background	The protein encoded by this gene is involved in cell motility. It is expressed in breast tissue and together with other proteins, it forms a complex with BRCA1 and BRCA2, thus is potentially associated with higher risk of breast cancer. Alternatively spliced transcript variants encoding different isoforms have been noted for this gene.
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

