



BDH2 Rabbit pAb

Catalog No	YP-Ab-17720
Isotype	IgG
Gene Name	BDH2
Alternative Names	DHRS6; EFA6R; SDR15C1; UCPA-OR; UNQ6308; PRO20933
Research Field	Cardiovascular
Product Categories	Primary antibody
Host	Rabbit
Reactivity	Human, Mouse
Application	WB, IHC-P
Dilution Ratio	WB: 1/500-1/1000 IHC: 1/50-1/100
Molecular Weight	Calculated MW: 27 kDa; Observed MW: 27 kDa
Clonality	Polyclonal Antibody
Immunogen	Full length fusion protein
Purification	Affinity Purified
Conjugation	Unconjugated
Modification	Unmodified
Form	Liquid
Buffer System	pH 7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Background	Dehydrogenase that mediates the formation of 2,5-dihydroxybenzoic acid (2,5-DHBA), a siderophore that shares structural similarities with bacterial enterobactin and associates with LCN2, thereby playing a key role in iron assimilation and homeostasis.
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

