



Mitofusin 1 Rabbit mAb

Catalog No	YP-Ab-17735
Isotype	IgG
Reactivity	Mouse,Rat,Hamster
Applications	WB,IP
Gene Name	Mfn1
Alternative Names	HR2; mKIAA4032; D3Ert265e; 2310002F04Rik; 6330416C07Rik
Research Field	Neuroscience
Product Categories	Primary antibody
Host	Rabbit
Molecular Weight	Calculated MW: 84 kDa; Observed MW: 84 kDa
Clonality	Monoclonal Antibody
Clonality No.	R05-3A8
Dilution	WB: 1/500-1/1000 IP: 1/20
Immunogen	A synthetic peptide of rat Mitofusin 1
Purification	Affinity Purified
Conjugation	Unconjugated
Modification	Unmodified
Form	Liquid
Buffer System	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Concentration	1 mg/ml
Purity	≥90%
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Background	Mitochondrial outer membrane GTPase that mediates mitochondrial clustering and fusion.
matters needing attention	Avoid repeated freezing and thawing!



Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

