



DLL3 Mouse mAb

Catalog No	YP-Ab-17736
Isotype	Mouse IgG2a
Reactivity	Human,Rat
Applications	WB,IHC-P,ICC/IF,FC,ELISA
Gene Name	DLL3
Alternative Names	SCDO1
Research Field	Neuroscience
Product Categories	Primary antibody
Host	Mouse
Molecular Weight	Calculated MW: 65 kDa;Observed MW: 65 kDa
Clonality	Monoclonal Antibody
Clonality No.	8Y7-D1-Y8
Dilution	WB: 1/500-1/1000 IHC: 1/100-1/200 IF: 1/50-1/200 FC: 1/50-1/100 ELISA: 1/10000
Immunogen	Purified recombinant fragment of human DLL3 expressed in E. Coli.
Purification	Affinity Purified
Conjugation	Unconjugated
Modification	Unmodified
Form	Liquid
Buffer System	Purified antibody in PBS with 0.05% sodium azide
Concentration	1 mg/ml
Purity	≥90%
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Background	This gene encodes a member of the delta protein ligand family. This family functions as Notch ligands that are characterized by a DSL domain, EGF repeats, and a transmembrane domain. Mutations in this gene cause autosomal recessive spondylocostal dysostosis 1. Two transcript variants encoding distinct isoforms have been identified for this gene.

**matters needing attention**

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

