



SELENOF Rabbit pAb

Catalog No	YP-Ab-17737
Isotype	IgG
Reactivity	Human,Mouse,Rat
Applications	WB
Gene Name	SELENOF
Alternative Names	selenoprotein F; SEP15
Research Field	-
Product Categories	Primary antibody
Host	Rabbit
Molecular Weight	18 kDa; 18kD(Calculated).
Clonality	Polyclonal Antibody
Clonality No.	-
Dilution	WB: 1/500-1/1000
Immunogen	Fusion protein of human SELENOF
Purification	Affinity Purified
Conjugation	Unconjugated
Modification	Unmodified
Form	Liquid
Buffer System	pH 7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Concentration	1 mg/ml
Purity	≥90%
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Background	May be involved in redox reactions associated with the formation of disulfide bonds (By similarity).
matters needing attention	Avoid repeated freezing and thawing!



Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

kDa

95—

72—

55—

36—

28—

17—

10—

