



TRMT10A rabbit pAb

□□	YP-Ab-17739
□□□	IgG
□□	WB
□□	Human;Mouse;Rat
□□	TRMT10A
□□□□	TRMT10A RG9MTD2
□□□	Synthesized peptide derived from human TRMT10A
□□□	This antibody detects endogenous levels of TRMT10A at Human, Mouse,Rat
□□	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
□□	Polyclonal, Rabbit,IgG
□□	WB 1:500-2000
□□□□	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
□□□	339 aa, 40 kDa
□□	S-adenosyl-L-methionine-dependent guanine N(1)-methyltransferase that catalyzes the formation of N(1)-methylguanine at position 9 (m1G9) in tRNAs . Probably not able to catalyze formation of N(1)-methyladenine at position 9 (m1A9) in tRNAs .
□□□□	Nucleus . Nucleus, nucleolus .
□□□□	Expressed in embryonic and fetal brain. It is expressed throughout the dorsal telencephalon at 8 and 11 weeks of gestation, with highest expression in ventricular zone and marginal zone. Detected in cerebellar cortex and nuclei, but not in dorsal telencephalon, at later stages.
□□	1 mg/ml
□□	-15°C to -25°C/1 year(Do not lower than -25°C)
□□□□□□	Avoid repeated freezing and thawing!
□□□□	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

