



PKG2 Rabbit pAb

Catalog No	YP-Ab-17741
Isotype	IgG
Reactivity	Human,Mouse,Rat
Applications	WB,IHC-P
Gene Name	PRKG2
Alternative Names	protein kinase; cGMP-dependent; type II; cGKII; PRKGR2
Product Categories	Primary antibody
Host	Rabbit
Molecular Weight	87kDa
Clonality	Polyclonal Antibody
Dilution	WB: 1/500-1/1000 IHC: 1/50-1/100
Immunogen	Synthetic peptide of human PRKG2
Purification	Affinity Purified
Conjugation	Unconjugated
Modification	Unmodified
Form	Liquid
Buffer System	pH 7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Concentration	1 mg/ml
Purity	≥90%
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Background	Crucial regulator of intestinal secretion and bone growth (By similarity).
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

