



MYH10 Rabbit mAb

Catalog No	YP-Ab-17742
Isotype	IgG
Reactivity	Human,Mouse
Applications	WB,IP
Gene Name	MYH10
Alternative Names	NMMHCB; NMMHC-IIB
Research Field	Signal Transduction
Product Categories	Primary antibody
Host	Rabbit
Molecular Weight	Calculated MW: 229 kDa; Observed MW: 229 kDa
Clonality	Monoclonal Antibody
Clonality No.	R05-6D8
Dilution	WB: 1/500-1/1000 IP: 1/20
Immunogen	A synthetic peptide of human non-muscle Myosin IIB/MYH10
Purification	Affinity Purified
Conjugation	Unconjugated
Modification	Unmodified
Form	Liquid
Buffer System	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Concentration	1 mg/ml
Purity	≥90%
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Background	Cellular myosin that appears to play a role in cytokinesis, cell shape, and specialized functions such as secretion and capping.
matters needing attention	Avoid repeated freezing and thawing!



Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

