



peptide YY rabbit pAb

□□	YP-Ab-17760
□□□	IgG
□□	PYY
□□□□	PYY
□□	Human;Mouse;Rat
□□	WB
□□□	Synthesized peptide derived from human peptide YY
□□□	This antibody detects endogenous levels of peptide YY at Human, Mouse,Rat
□□	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
□□	Polyclonal, Rabbit,IgG
□□	WB 1:500-2000
□□□□	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
□□□	97 aa, 11 kDa
□□	This gut peptide inhibits exocrine pancreatic secretion, has a vasoconstrictory action and inhibits jejunal and colonic mobility.
□□□□	Secreted.
□□	1 mg/ml
□□	-15°C to -25°C/1 year(Do not lower than -25°C)
□□□□□□	Avoid repeated freezing and thawing!
□□□□	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images