



Plexin A1 Rabbit mAb

Catalog No	YP-Ab-17768
Isotype	IgG
Reactivity	Human,MouseRat
Applications	WB,ICC/IF,FC
Gene Name	PLXNA1
Alternative Names	NOV; NOV; Plexin-A1; Plexin1; PlexinA1; PLXN1; Plxna1
Research Field	Neuroscience
Product Categories	Primary antibody
Host	Rabbit
Molecular Weight	Calculated MW: 211 kDa; Observed MW: 211 kDa
Clonality	Monoclonal Antibody
Clonality No.	R08-2J4
Dilution	WB: 1/500-1/1000 IF: 1/50-1/200 FC: 1/50-1/100
Immunogen	A synthesized peptide derived from human Plexin A1
Purification	Affinity Chromatography
Conjugation	Unconjugated
Modification	Unmodified
Form	Liquid
Buffer System	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Concentration	1 mg/ml
Purity	≥90%
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Background	Plays a role in axon guidance, invasive growth and cell migration.
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images