



ApoA2 Rabbit mAb

Catalog No	YP-Ab-17769
Isotype	IgG
Reactivity	Human,Rat
Applications	WB,ICC/IF,FC
Gene Name	APOA2
Alternative Names	APO A2; Apo AII; APOA 2; ApoA II; APOA2; APOAII; Apolipoprotein A II
Research Field	Cardiovascular
Product Categories	Primary antibody
Host	Rabbit
Molecular Weight	Calculated MW: 11 kDa; Observed MW: 11 kDa
Clonality	Monoclonal Antibody
Clonality No.	R09-8W8
Dilution	WB: 1/500-1/1000 IF: 1/50-1/200 FC: 1/50-1/100
Immunogen	A synthesized peptide derived from human Apolipoprotein A II
Purification	Affinity Chromatography
Conjugation	Unconjugated
Modification	Unmodified
Form	Liquid
Buffer System	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Concentration	1 mg/ml
Purity	≥90%
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Background	May stabilize HDL (high density lipoprotein) structure by its association with lipids, and affect the HDL metabolism.
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images