



SLC2A9 rabbit pAb

□□	YP-Ab-17784
□□	WB
□□	Human
□□□	IgG
□□	SLC2A9
□□□□	SLC2A9 GLUT9
□□□	Synthesized peptide derived from human SLC2A9
□□□	This antibody detects endogenous levels of SLC2A9 at Human
□□	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
□□	Polyclonal, Rabbit,IgG
□□	WB 1:500-2000
□□□□	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
□□□	540 aa, 59 kDa
□□	Urate transporter, which may play a role in the urate reabsorption by proximal tubules . Does not transport glucose, fructose or galactose .
□□□□	[Isoform 1]: Basolateral cell membrane ; Multi-pass membrane protein .; [Isoform 2]: Apical cell membrane ; Multi-pass membrane protein. Basolateral cell membrane ; Multi-pass membrane protein .
□□□□	Most strongly expressed in basolateral membranes of proximal renal tubular cells, liver and placenta. Also detected in lung, blood leukocytes, heart skeletal muscle and chondrocytes from articular cartilage. Isoform 2 is only detected in the apical membranes of polarized renal tubular cells and placenta. Isoform 1 and isoform 2 are detected in kidney membrane (at protein level).
□□	1 mg/ml
□□	-15°C to -25°C/1 year(Do not lower than -25°C)
□□□□□□	Avoid repeated freezing and thawing!
□□□□	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.