



SOX10 Rabbit mAb

Catalog No	YP-Ab-17800
Isotype	IgG
Reactivity	Human,Rat
Applications	WB
Gene Name	SOX10
Alternative Names	Transcription factor SOX-10
Research Field	Neuroscience
Product Categories	Primary antibody
Host	Rabbit
Molecular Weight	Calculated MW: 50 kDa; Observed MW: 62 kDa
Clonality	Monoclonal Antibody
Clonality No.	R09-7D8
Dilution	WB: 1/500-1/1000
Immunogen	A synthetic peptide of human SOX10
Purification	Affinity Purified
Conjugation	Unconjugated
Modification	Unmodified
Form	Liquid
Buffer System	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Concentration	1 mg/ml
Purity	≥90%
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Background	Sox10 is a member of the SOX (SRY-related HMG-box) family of transcription factors involved in the regulation of embryonic development and the determination of cell fate. Sox10 is an important regulator of neural crest and peripheral nervous system development.



matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

