



Bcl2 (1D8) Mouse mAb

Catalog No	YP-Ab-17814
Isotype	IgG1
Reactivity	Human,Rat,Mouse
Applications	WB
Gene Name	BCL2
Alternative Names	BCL2; Apoptosis regulator Bcl-2
Research Field	Cell Biology
Product Categories	Primary antibody
Host	Mouse
Molecular Weight	Calculated MW: 26 kDa; Observed MW: 26 kDa
Clonality	Monoclonal Antibody
Clonality No.	1D8-1G7-1G3
Dilution	WB: 1/500-1/1000
Immunogen	Synthetic peptide conjugated to KLH.
Purification	Affinity Purified
Conjugation	Unconjugated
Modification	Unmodified
Form	Liquid
Buffer System	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
Concentration	1 mg/ml
Purity	≥90%
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Background	This gene encodes an integral outer mitochondrial membrane protein that blocks the apoptotic death of some cells such as lymphocytes. Constitutive expression of BCL2, such as in the case of translocation of BCL2 to Ig heavy chain locus, is thought to be the cause of follicular lymphoma.



matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

