

Mouse Anti Human sFLT-1 mAb-Clone 2

Catalog No□

Quantity: Specified on product label

Buffer: PBS,PH7.4

Storage: Store at 2-8C short term and -80°C long term.Storage

It is recommended to aliquot needed volumes and store at -80C

Avoid repeated freeze/thaw cycles

QC Results

Characteristics	Method	Requirements	Results
Appearance	Visual observation	Clear Liquid	Qualified
Concentration	A280 with the extinctioncoefficient 1.43	>2.00mg/mL	8.26mg/mL
Purity	SDS-PAGE	295%	>95%
Calculated Molecular Weight			