



ACSF2 Rabbit pAb

Catalog No	YP-Ab-17824
Isotype	IgG
Reactivity	Human,Mouse,Rat
Applications	WB,IHC-P
Gene Name	ACSF2
Alternative Names	acyl-CoA synthetase family member 2; ACSMW; AVYV493
Research Field	-
Product Categories	Primary antibody
Host	Rabbit
Molecular Weight	615 aa, 68 kDa
Clonality	Polyclonal Antibody
Clonality No.	-
Dilution	WB: 1/500-1/1000 IHC: 1/50-1/100
Immunogen	Fusion protein of human ACSF2
Purification	Affinity Purified
Conjugation	Unconjugated
Modification	Unmodified
Form	Liquid
Buffer System	pH 7.4 PBS, 0.05% NaN3, 40% Glycerol
Concentration	1 mg/ml
Purity	≥90%
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Background	Acyl-CoA synthases catalyze the initial reaction in fatty acid metabolism, by forming a thioester with CoA.
matters needing attention	Avoid repeated freezing and thawing!



Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

