



FOXO6 rabbit pAb

□□	YP-Ab-17835
□□□	IgG
□□	WB
□□	Human;Mouse
□□	FOXO6
□□□□	FOXO6
□□□	Synthesized peptide derived from human FOXO6
□□□	This antibody detects endogenous levels of FOXO6 at Human, Mouse
□□	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
□□	Polyclonal, Rabbit,IgG
□□	WB 1:500-2000
□□□□	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
□□□	51 kDa
□□	Transcriptional activator.
□□□□	Cytoplasm . Nucleus . When phosphorylated, translocated from nucleus to cytoplasm. High nuclear localization after stimulation with growth factors (By similarity). .
□□	1 mg/ml
□□	-15°C to -25°C/1 year(Do not lower than -25°C)
□□□□□□	Avoid repeated freezing and thawing!
□□□□	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images