



PI 3 Kinase p150 Rabbit mAb

Catalog No	YP-Ab-17849
Isotype	IgG
Reactivity	Human,Mouse,Rat
Applications	WB,IHC-P,ICC/IF,FC
Gene Name	PIK3R4
Alternative Names	P150; Phosphoinositide 3 kinase regulatory subunit 4; PI3 kinase p150 subunit; PIK3R4; VPS15
Research Field	Signal Transduction
Product Categories	Primary antibody
Host	Rabbit
Molecular Weight	Calculated MW: 153 kDa; Observed MW: 150 kDa
Clonality	Monoclonal Antibody
Clonality No.	R07-6N5
Dilution	WB: 1/500-1/1000 IHC: 1/50-1/100 IF: 1/50-1/200 FC: 1/50-1/100
Immunogen	A synthesized peptide derived from human PI 3 Kinase R4
Purification	Affinity Purified
Conjugation	Unconjugated
Modification	Unmodified
Form	Liquid
Buffer System	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Concentration	1 mg/ml
Purity	≥90%
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Background	Regulatory subunit of the PI3K complex. May regulate membrane trafficking late in the endocytic pathway.
matters needing attention	Avoid repeated freezing and thawing!



Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

