



UBE2C rabbit pAb

□□	YP-Ab-17997
□□□	IgG
□□	WB
□□	Human;Mouse
□□	UBE2C
□□□□	UBE2C UBCH10
□□□	Synthesized peptide derived from human UBE2C
□□□	This antibody detects endogenous levels of UBE2C at Human, Mouse
□□	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
□□	Polyclonal, Rabbit,IgG
□□	WB 1:500-2000
□□□□	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
□□□	179 aa, 20 kDa
□□	Accepts ubiquitin from the E1 complex and catalyzes its covalent attachment to other proteins. In vitro catalyzes 'Lys-11'- and 'Lys-48'-linked polyubiquitination. Acts as an essential factor of the anaphase promoting complex/cyclosome (APC/C), a cell cycle-regulated ubiquitin ligase that controls progression through mitosis. Acts by initiating 'Lys-11'-linked polyubiquitin chains on APC/C substrates, leading to the degradation of APC/C substrates by the proteasome and promoting mitotic exit.
□□	1 mg/ml
□□	-15°C to -25°C/1 year(Do not lower than -25°C)
□□□□□□	Avoid repeated freezing and thawing!
□□□□	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

