



SLC25A10 rabbit pAb

□□	YP-Ab-18081
□□□	IgG
□□	WB
□□	Human;Mouse
□□	SLC25A10
□□□□	SLC25A10 DIC
□□□	Synthesized peptide derived from human SLC25A10
□□□	This antibody detects endogenous levels of SLC25A10 at Human, Mouse
□□	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
□□	Polyclonal, Rabbit,IgG
□□	WB 1:500-2000
□□□□	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
□□□	296 aa, 32 kDa
□□	Involved in translocation of malonate, malate and succinate in exchange for phosphate, sulfate, sulfite or thiosulfate across mitochondrial inner membrane.
□□□□	Mitochondrion inner membrane; Multi-pass membrane protein.
□□□□	Present in high amounts in liver and kidney, and at lower levels in all the other tissues analyzed.
□□	1 mg/ml
□□	-15°C to -25°C/1 year(Do not lower than -25°C)
□□□□□□	Avoid repeated freezing and thawing!
□□□□	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images